



Order: (888)-282-5810 (Phone)  
(818)-707-0392 (Fax)  
[order@abiocode.com](mailto:order@abiocode.com)  
Web: [www.Abiocode.com](http://www.Abiocode.com)

---

## Protein Extracts for GLR3.2 Analysis

**Cat#: L3326**

**Lot#: Refer to vial**

**Application: WB**

**Quantity: 50 ul WT and 50 ul *GLR3.2* mutant extracts**

### **Product Description:**

Protein extracts were prepared from wild type (WT, Columbia) or *GLR3.2*-knockdown mutant (Stock#: SALK\_125952C, intergenic T-DNA insertion) leaves. Plant tissues were ground to fine powders and mixed with 4 x SDS sample buffer. This product can be used as controls for Western blot analysis of *GLR3.2* (At4g35290) protein (Uniprot ID: Q93YT1).

### **Storage Condition:**

Keep at -20°C. Stable for 6 months from date of receipt.

### **Usage Instruction:**

For Western blot analysis, incubate protein extracts at 90°C for 5 min and load 10 ul of each extract per lane.